

Standard operating procedure (SOP) for human *ex vivo* skin microdialysis

A protocol from the EAACI Task Force on Skin Microdialysis

Aim

The aim of this document, which is a product of the EAACI Task Force on Skin Microdialysis, is to describe a standardized procedure for sampling of soluble molecules from human skin *ex vivo* using microdialysis.

Skin microdialysis (SMD) facilitates recovery of endogenous and exogenous molecules from the extracellular compartment and has been used, primarily *in vivo*, to study many dermal processes, including skin diseases (e.g. atopic dermatitis and psoriasis), percutaneous absorption, hypersensitivity reactions and cutaneous inflammation.

Please refer to the corresponding position paper from the task force for a more detailed description of SMD and its applications. The position paper also features pictures of a SMD experiment ([Fig. 2](#)).

Timeline

The duration of a SMD procedure is highly dependent on the scope of the study and its setup. However, time estimates are listed for each part of the experiment below:

- | | |
|--|--|
| A. Preparation of collection vials ¹ | <i>Allow for 2-3 hours depending on the number of collection vials</i> |
| B. Preparation of pump and microdialysis probes | <i>Allow for 1-2 hours</i> |
| C. Preparation of human <i>ex vivo</i> skin and sampling | |
| - Skin preparation | <i>Allow for 30-60 minutes</i> |
| - Probe insertion | <i>Allow for 15-30 minutes for insertion per probe</i> |
| - Microdialysis sampling | <i>Minutes to days depending on the setup</i> |
| D. Post-sampling procedure | <i>Allow for 15 minutes</i> |

Note

Make sure that you are acquainted with the equipment required before commencing SMD experiments.

¹ Preparation of collection vials can be performed the day before the actual SMD sampling is conducted or during priming of the microdialysis probes (depending on the number of probes and vials).

Materials (must be available when commencing the experiment)

Materials	Use and notes
Catheter connectors	To connect the probe inlet tubing to syringes mounted in the pump
Collection vials	For collection of dialysates during sampling, e.g. PCR tubes or glass capillaries ²
Forceps	For handling of probes
Gloves	For aseptic handling
Guide cannulas	For intradermal insertion of microdialysis probes, e.g. 21G/23G syringe needles depending on the probe outer diameter or 25G needles with attached (crimped) probes
Ice	<i>Optional:</i> For stabilization of samples during sampling/prior to storage
Marker	To mark the treatment area and the probe insertion sites on the skin
Microdialysis probes	Self-made or commercially available; linear or concentric. Low molecular weight cut-off (MWCO) (e.g. 2-5 kDa) for recovery of small molecules and high MWCO (e.g. 100-3000 kDa) probes for recovery of larger molecules
Microinfusion pump	For perfusion of microdialysis probes, with stable flow rates between 0.1-5 μ l/min
Paper towel	Wetted with saline to keep the <i>ex vivo</i> skin specimens hydrated during the entire procedure
Parafilm	To cover collection vials in order to minimize evaporation during sampling
Perfusate	For perfusion (often a sterile physiological solution, but should be adapted to fit the target molecules ³)
Ruler	To measure the intradermal insertion length of microdialysis probes
Saline (0.9% NaCl solution)	For flushing of guide cannulas to remove dead cells and blood and to ease probe insertion, to rinse <i>ex vivo</i> skin specimens, and to maintain a moist environment

² The collection vials must be appropriate for the study undertaken and should take into account the characteristics of the analyte. Important properties to consider when choosing the type of collection vial: Material (e.g. low adherence plastic tubes, possibly containing a stabilizing agent to diminish analyte breakdown, or glass capillaries), volume, label-ability, storage/freeze-ability, type of lid/seal, compatibility with refrigerated fraction collectors, etc.

³ The perfusate composition should be modified based on the characteristics of the target analyte, e.g. if highly lipophilic, a lipid emulsion can be used. Additionally, when using high MWCO probes to recover large molecules a colloid may be added to counteract fluid leakage from the probe caused by the larger pores in the microdialysis membrane.

Scale	To weigh collection vials
Scalpel	To excise <i>ex vivo</i> skin specimens from the parent tissue sample
Scissors	To adjust the probe outlet
Styrofoam	To pin <i>ex vivo</i> skin specimens onto for better handling
Surgical scissors	To remove subcutaneous fat from <i>ex vivo</i> skin specimens
Surgical tape	To fixate microdialysis probes during priming and sampling
Syringes (flush)	To flush probes and guide cannulas
Syringes (pump)	To contain the perfusate (gas-tight or disposable plastic syringes compatible with the pump and volume delivery required)
Ultrasound scanner	<i>Optional:</i> To assess the intradermal probe depth

Protocol

A) Preparation of collection vials⁴ (allow for 2-3 hours depending on the number of vials)

Step	Procedure
1	Choose appropriate collection vials ²
2	Label collection vials
3	Weigh individual collection vials and note the weight before sampling ^{5,6}
4	Adapt the outlet tubing to the collection vial (e.g. using a plastic stopper for glass capillaries or cover collection vials with Parafilm and punch a small hole for probe insertion using a needle or a lancet)

B) Preparation of pump and microdialysis probes (allow for 1-2 hours)

Step	Procedure
1	If necessary, adjust the length of the (linear) microdialysis membrane ⁷
2	Prepare the pump syringes for perfusion by filling them with perfusate ⁸ (e.g. Ringer-lactate + 1% human serum albumin (HSA) ⁹)
3	Mount the syringes in the microinfusion pump ¹⁰
4	Gently flush the probes with perfusate (with a sufficient volume to ensure that the probe only contains perfusate, e.g. 200 µl) using a flush button on the pump or manually while following the fluid through the probe ¹¹
5	Confirm flow through the probe and ensure that it does not leak

⁴ Preparation of collection vials can be performed the day before the actual SMD sampling is conducted or during priming of the microdialysis probes (depending on the number of probes and vials).

⁵ The fluid recovery can be monitored by weighing collection vials before and after sampling. This is especially important when using probes with a high MWCO.

⁶ Weighing and covering of collection vials may be omitted if sampling is carried out at higher flow rates and for shorter time intervals, as this decreases the risk of the sample evaporating.

⁷ Be careful not to damage the delicate membrane when adjusting its length (make sure to use sharp scissors).

⁸ Avoid formation of bubbles in probes and syringes, especially when the perfusate contains proteins (e.g. HSA).

⁹ When using high MWCO probes colloid additives should be included in the perfusate. For instance, HSA-enriched perfusates can be used for sampling of larger molecules such as cytokines, as it minimizes non-specific adsorption to probe components and acts as a carrier protein.

¹⁰ Only calibrated microinfusion pumps should be used (the flow rates used must be validated).

¹¹ Be careful not to push out the central guide wire when flushing linear probes.

6	Connect the probe inlet tubing to the syringes mounted in the pump using catheter connectors ¹²
7	Prime microdialysis probes (for up to 1 hour if the perfusate contains 1% HSA to minimize non-specific adsorption) ¹³ (see Fig. 2A in the position paper)

C) Preparation of human *ex vivo* skin and sampling¹⁴ (allow for 30-60 minutes of skin preparation and 15-30 minutes for insertion of each probe – the actual sampling may take from minutes to days)

Step	Procedure
1	Pin the parent <i>ex vivo</i> tissue sample onto Styrofoam to allow for visual inspection and handling of the skin
2	Excise appropriate skin specimens from the parent tissue sample using a scalpel ¹⁵
3	Remove any subcutaneous fat (e.g. using a scalpel and/or a scissor) without causing unnecessary trauma to the tissue
4	Pin the excised skin specimens (dermal side down) onto Styrofoam with a wet paper towel in between to maintain a moist environment
5	Mark the probe entry and exit points (for linear probes) on the <i>ex vivo</i> skin specimens ¹⁶ (e.g. placing the probes 2 cm apart ¹⁷)
6	Insert the guide cannulas (either 21G or 23G syringe needles depending on the outer probe diameter) spanning the planned intradermal length (e.g. 2 cm) ^{18,19}
7	Flush each guide cannula inside the skin with sterile saline to ease insertion of the probe
8	Insert the primed microdialysis probes through the tip of the guide cannulas while the pump is running. Repeat for all probes (see Fig. 2D in the position paper)
9	Ensure that dialysate appears from the outlet tubing and gently pull out the guide cannulas, thus leaving the probe in place

¹² Be careful not to introduce bubbles in any part of the system.

¹³ High MWCO probes may be tilted downwards (approximately 45°) during priming and sampling to ensure a proper flow.

¹⁴ Gloves should be worn during the entire procedure to protect the operator and the microdialysis membranes.

¹⁵ Avoid areas with visible local bleeding, lesions, scars or stretch marks.

¹⁶ Make sure there are no scars, lesions or irregularities at the sites of probe insertion.

¹⁷ In general, probes should be placed at least 1 cm apart to avoid spillover, but the interprobe distance must be adapted to the analyte diffusibility in the tissue. Furthermore, probes should be distanced at least 1 cm from the edge of the skin specimen to avoid spillover from the excision trauma.

¹⁸ The insertion depth of the microdialysis probes must be kept as consistent as possible, for which reason the same operator should insert all guide cannulas throughout a study. The probe depth can be assessed using an ultrasound scanner.

¹⁹ The guide cannulas must be inserted evenly without penetrating the skin barrier apart from at the point of probe entry and at the exit point, if linear probes are used.

10	Gently fix the probe onto the skin using a piece of surgical tape placed on the inlet tubing immediately proximal to the membrane (see Fig. 2G in the position paper)
11	Tilt the Styrofoam with <i>ex vivo</i> skin specimens pinned onto it to approximately 45° to ensure optimal flow ²⁰
12	Place the collection vials at the outlet of the microdialysis probes and secure the vials using pieces of surgical tape if needed ²¹ (see Fig. 2I in the position paper)
13	Start SMD sampling, e.g. after performing relevant experimental procedures such as injection of test substances ^{22,23}
14	When the SMD sampling is completed, reweigh collection vials
15	Store the dialysates according to the analyte of interest (e.g. at –80 °C for subsequent protein analysis)

D) Post-sampling procedure (allow for 5-15 minutes)

Step	Procedure
1	Discard waste materials (such as needles and tissue samples) in accordance with local procedures

Abbreviations

HSA	Human serum albumin
MWCO	Molecular weight cut-off
SMD	Skin microdialysis

²⁰ High MWCO probes must be tilted downwards during priming and sampling to ensure a proper flow. This can be omitted for low MWCO probes.

²¹ If labile analytes are studied, the collection vials should be cooled during or immediately after sampling.

²² Equilibration time may be added prior to sampling to reduce the impact of trauma from probe insertion.

²³ Remember to accommodate for sampling dead space in the collection protocol.