

# BODY COMPOSITION ANALYZER

# **SC-331S**

## Instruction manual

SC-331S

SC-331S with Column Kit

### <Usage Conditions>

Temperature Range for Use : 0°C – 35°C  
Relative Humidity : 30% – 80%  
(without condensation)  
Max Altitude : 2,000m ASL  
Pressure Range of Environment : 86kPa - 106kPa

### <Storage Conditions>

Temperature Range of Environment : -10°C – 60°C  
Range of Relative Humidity : 10% – 90% (without condensation)  
To avoid malfunctions, avoid storing the equipment where there is direct sunlight, significant temperature changes, the risk of dampness, a large amount of dust, in the vicinity of fires, or where there is the risk of receiving vibrations or shocks.

Please read this Instruction Manual carefully and keep it handy for future reference.

## Intended use

- This equipment can be used in the screening of certain adult diseases and conditions related to body weight and composition.
- It can be used in the monitoring and prevention of conditions caused by excessive deposits of fat tissue such as diabetes, hyperlipidemia, cholelithiasis and fatty liver.
- It can be used in the monitoring of changes in individuals' body composition related to differences in the ratio of fat tissue to lean.
- It can be used to assess the effectiveness of individuals' nutrition and exercise programmers, both for health and physical fitness.
- The TANITA Body Composition Analyzer is indicated for use in the measurement of weight and impedance, and the estimation of body mass index (BMI), total body fat percent, total body water percent and weight, muscle mass (skeletal and smooth), physique rating, bone mass, visceral fat rating with healthy range, basal metabolic rate (BMR), metabolic age, and target body fat percent with predicted weight and fat mass, using BIA (Bioelectrical Impedance Analysis).  
The device is indicated for use for healthy children 5-17 years old and healthy adults with active, moderately active, to inactive lifestyles.

## Efficacy

1. This product is simple to use, and requires no specialized facilities or expertise to take measurements.
2. Measurements can be taken quickly and easily, causing minimal inconvenience to the patient during measurements.

## Contents

### Before use

**Safety Notes** .....4

**Product Assembly and Components**.....8

**Preparation**.....11

**Various settings**.....12

### How to use

**Operating Instructions**.....28

- When using as a body composition analyzer .....28
- Target body fat % .....38
- when using as a scale .....39

**Various criteria**.....42

### When necessary

**Troubleshooting** .....48

**Connection with a personal computer** .....50

**Technical Notes** .....54

**The New Regression Formula for Basal Metabolic Rate (BMR)**.....58

**Specifications** .....60

# For Your Safety

The following outlines precautionary measures to be taken to avoid injury to the users of this device and others, and to prevent damage to property. Please familiarise yourself with the contents to ensure the safe operation of this equipment.

**Warning** Ignoring instructions highlighted with this mark could result in death or severe injury.

**Caution** Ignoring instructions highlighted with this mark could result in injury or damage to property.

This mark indicates actions that are prohibited.

This mark indicates instructions that must always be adhered to.

## Warning

**This device must not be used on subjects with pacemakers or other mechanical implants.**

This device passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious repercussions.

**Do not handle the plug with a wet hand.**

Such action may result in electric shock, fire, or leakage.

**Keep this equipment away from flammable gas and oxygen rich environment.**


**No modification of this equipment is allowed.**

There is a risk of electrocution or injury. Nor can precise analysis be guaranteed.

**Do not use double adapters.**

May result in fire.

**Use only a correctly wired outlet.**

Only genuine cables and equipments are allowed.

Always...

## Caution

**Keep off wet circumstance.**

**Avoid using on subjects with allergies to metals.**

Allergic reactions may be caused by the stainless steel used in the electrodes of this device.

**Do not jump on the equipment.**

**Do not lean on the equipment.**

Do not use this equipment near other products that emit electromagnetic waves.

**Do not insert fingers into gaps and holes.**

**Avoid impacts on the display.**

The screen panel may split and cause injury.

**Ensure disabled persons do not attempt to use the equipment by themselves.**

A companion must perform the measurements.

Prohibited

Always...

**Be sure to clean the scale platform with appropriate disinfectant after each use.**

**For your safety to avoid the risk of electric shock and to secure the accuracy, keep a clearance with patients during measuring.**

**Continually monitor both the subject and the equipment for anomalies.**

Should an anomaly with the subject or the equipment be discovered, take appropriate action, such as stopping the equipment, while ensuring the safety of the subject.

**Unplug the equipment from the wall outlet when it will not be in use for long periods of time.**

**Do not interpret analysis results yourself (including evaluating measurements and formulating exercise programmes based on results, etc.).**

Weight loss and exercise based on self-analysis could be detrimental to you health. Always follow the advice of a qualified professional.

**This equipment is designated a Class B IT equipment (mainly for systems intended to be used in domestic environments) and CE (EMC), but weak equipment may be affected.**

If you connect your computer or peripherals, please use the compliance with IEC60601-1 (EN60601-1). Or if using the compliance with IEC60950 (EN60950), power must be supplied from the medical isolating transformer. Keep at a distance of 1.5m between each equipment during operation. If you ignore this requirement, may cause electrical shock of patients or operation.

## For Accurate Measurements

### Avoid measuring after hard exercise.

May result in measurement errors. Please take measurements after sufficient rest.

### Avoid measuring after over-eating or over-drinking, and when severely dehydrated.

May result in measurement errors. For greater accuracy, avoid using directly after waking up. Use at the same time on each occasion, at least three hours after eating.

### Do not take measurements while using transmitters, such as mobile phones, which may affect readings.

### Use the device under the same conditions and in the same position as much as possible to track changes.

Readings are greatly affected by the level of hydration and position of the body. Please use at the same time each day, under the same conditions and in the same body position.

### Avoid measuring in locations of greatly differing temperature.

May result in measurement errors. When the equipment has been moved to a location  $\pm 20^{\circ}\text{C}$ , allow to stand for at least 2 hours before using.

### Bare feet should be placed correctly on the electrode panel for measurement.

Also, make sure the soles of feet are free of excess dirt, as this may also act as a barrier to the mild current.

### Use in a stable location.

Errors in measurements may occur when the device is used in an unstable location.


# Product Assembly and Components

Eng


Before use  
(Product Assembly and Components)

## SC-331S

### Control Box


### Platform


### Accessories

- Instruction manual (This manual)
- Quick reference guide
- CD-ROM (USB driver)
- AC adapter
- Power cord

- Dropper (1 item)


- Printer paper

(Thermal paper, roll diameter:  
2.2in / 55mm, roll length:  
approx. 112ft / 34m)

\* Please contact the agent from  
which you have purchased  
the product for details.


## SC-331S with Column Kit


### Accessories

- Instruction manual (This manual)
- Quick reference guide
- CD-ROM (USB driver)
- Assembling guide
- AC adapter
- Power cord

- Hexagonal wrench (1 item)
- Hexagon socket head bolts  
(M5L12) (4 items)

- Dropper (1 item)

- Printer paper

(Thermal paper, roll diameter:  
2.2in / 55mm, roll length:  
approx. 112ft / 34m)

\* Please contact the agent from  
which you have purchased  
the product for details.



**Caution** Make sure you place the Platform on a stable, level surface. If the Platform is not stable because not all the feet are on the surface, for example, there is a risk of stumbling or inaccurate measurement.

**Caution** Make sure you place the Platform on a stable, level surface. If the Platform is not stable because not all the feet are on the surface, for example, there is a risk of stumbling or inaccurate measurement.

Eng

Before use  
(Product Assembly and Components)

### Warning

- To avoid electric shock, do not insert or remove the plug with wet hands.
- To avoid electric shock, do not use the equipment near water.

### Caution

- To avoid measurement error, do not measure while using equipment that generates radio waves, such as cellular phone.
- Use only the original AC adapter. Using AC adapters other than the original one may cause malfunction, smoke or fire.

### Caution

- Please change the printer paper when red lines appear along the sides of the paper.
- Please pay careful attention to avoid injury from the sharp edge.
- Please turn off the machine before clearing Paper jams.


### Setting of the printer paper

**1** Press **ON/OFF** to turn on the power.

- After all lamps light up, the model number is displayed as **331**, and **00 kg** is displayed.

**2** (1) Press **OPEN**.  
(2) Remove the printer cover.


**3** (3) Set up the printer paper

- Remove the adhesive of the printer paper and draw it out approx. 4in / 10cm.


**4** Return the printer cover.

- In the case that **COPE n** is displayed, ⇨ The printer cover is open, so please close it again properly (☞ page 49).

**5** Press the **FEED**, and cut off the excessive paper.

- In the case that the automatic cutting is set to "OFF," the automatic cutting does not function (☞ page 16).

**6** Setting completed.

# Various settings

## Setting methods

Eng

Before use  
(Various settings)


Eng

Before use  
(Various settings)

### Call up the setting item selection screen.

#### 1 Press **ON/OFF** to turn on the power.

- After all lamps light up, the model number is displayed, and **00 kg** is displayed.


#### 2 Press **Mode Setting**

- The setting item input screen is displayed.

#### Note

- 39 – 64 are the setting items related to print out items (page 24).
- When the various settings are all completed,
  - ⇒ Press **Mode Setting** on the “Setting item selection” screen (it returns to the preset tare input screen).

#### Setting items

1	Date and time (page 14)
2	Number of sheets to print Body composition analyzer (page 16)
3	Number of sheets to print Weight only mode (page 15)
4	Auto Cutting of the printer paper (page 16)
5	Beep sound (page 16)
6	Display Fat % Healthy range (page 17)
7	ID Number (page 17)
8	Measurement flow (page 18)
9	Athletic mode (page 18)
10*	Height increment (page 19)
11	Automatic determination time (page 19)
17	Wrestler mode (page 20)
18	Target body fat % (page 21)
19	Select language (page 21)
20	Print item preset (page 22)

\* Invalid in “lb” mode.

#### Print item setting ON/OFF (page 24)

39	TANITA Logo
40	Category name
41	Date and Time
42	Serial number
43	Memo space
44	ID Number
45	Fat mass
46	Fat free mass
47	Muscle mass
48	Total body water
49	Total body water %
50	Bone mass
51	Basal metabolic rate
52	Metabolic age
53	Visceral fat rating
54	BMI
56	Ideal body weight
57	Degree of obesity
58	Desirable range Body fat
59	Indicator - Fat %
60	Indicator - BMI
61	Indicator - Visceral fat rating
62	Indicator - Muscle mass
63	Indicator - BMR
64	Physique rating

#### Note

- \* When various settings are continuously carried out,
  - ⇒ press each number to set.
- \* The set contents are memorized until they are changed next time.


\* When various settings are all completed,

⇒ press **Mode Setting** on the “setting item selection screen” (it returns to the preset tare input screen).

## Register the date and time (Continued from page 12).

**3** Press **1** and then press **Enter / Next**.

- The "date and time" setting screen is displayed.


**4** Input the date and time

- Input the year, month, day, hour and minute in order with 2 digits.

Example: In the case to input: 9:47 am January 21<sup>st</sup>, 2009


0 9 0 1 2 1 0 9 4 7

- Note**
- To input a number with 1 digit (0 - 9), first press **0**.
  - To correct the input value, or cancel the input, press **CE** (the input is deleted).
  - To end inputting in midstream, press **Enter / Next**.

**5** After inputting all the items, press **Enter / Next**.

- It returns to the "setting item selection" screen.


**Note**

- \* When various settings are continuously carried out, press each number to set.
- \* The set contents are memorized until they are changed next time.

## Set the number of sheets to print for the body composition Analyzer (Continued from page 12).



**3** Press **2** and then press **Enter / Next**.

- The "number of sheets to print for the body composition analyzer" setting screen is displayed.

**4** Input the number of sheets to print.

- Note**
- The default is "1" (input range: 0 - 3).
  - In the case that the number of sheets to print is set to "0" for both the body composition analyzer and the scale, **FEED** key does not function.
  - To correct the input value, or cancel the input, press **CE** (the input is deleted).


**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.

## Set the number of sheets to print for the scale (Continued from page 12).



**3** Press **3** and then press **Enter / Next**.

- The "number of sheets to print for the scale" setting screen is displayed.

**4** Input the number of sheets to print.

- Note**
- The default is "1" (input range: 0 - 3).
  - In the case that the number of sheets to print is set to "0" for both the body composition analyzer and the scale, **FEED** key does not function.
  - To correct the input value, or cancel the input, press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.


- It returns to the "setting item selection" screen.

- \* When various settings are all completed, press **Mode Setting** on the "setting item selection screen" (it returns to the preset tare input screen).


## Set ON / OFF of the automatic cutting of the printer paper (Continued from page 12).


- 3** Press **4** and then press **Enter / Next**.
- The "ON or OFF of automatic cutting" selection screen is displayed.

- 4** Set ON or OFF of the automatic cutting.
- Note**
- The default is "0.off." ("1. on" for valid, "0. off" for invalid.)
  - To correct the input value, or cancel the input, ⇨ press **CE** (the input is deleted).


- 5** After inputting the numeric value, press **Enter / Next**.
- It returns to the "setting item selection" screen.

## Set ON / OFF of the beep sound (Continued from page 12).


- 3** Press **5** and then press **Enter / Next**.
- The "ON or OFF of the beep sound" selection screen is displayed.

- 4** Set ON or OFF of the beep sound.
- Note**
- The default is "1.on." ("1. on" for valid, "0. off" for invalid.)
  - To correct the input value, or cancel the input, ⇨ press **CE** (the input is deleted).

- 5** After inputting the numeric value, press **Enter / Next**.
- It returns to the "setting item selection" screen.


**Note**


- \* When various settings are continuously carried out, ⇨ press each number to set.
- \* The set contents are memorized until they are changed next time.

## Set ON / OFF of the Fat % Healthy range display (Continued from page 12).

- 3** Press **6** and then press **Enter / Next**.
- The "ON or OFF of the Fat % Healthy range" selection screen is displayed.


- 4** Set ON or OFF of the Fat % Healthy range display when using the body composition


- Note**
- The default is "1.on." ("1. on" for valid, "0. off" for invalid.)
  - To correct the input value, or cancel the input, ⇨ press **CE** (the input is deleted).

- 5** After inputting the numeric value, press **Enter / Next**.
- It returns to the "setting item selection" screen.


## Set with or without an ID (Continued from page 12).


- 3** Press **7** and then press **Enter / Next**.
- The "with or without an ID" setting screen is displayed.

- 4** Set with or without an ID.
- Note**
- The default is "0.off." ("1. on" for valid, "0. off" for invalid.)
  - To correct the input value, or cancel the input, ⇨ press **CE** (the input is deleted).

- 5** After inputting the numeric value, press **Enter / Next**.
- It returns to the "setting item selection" screen.

- \* When various settings are all completed, ⇨ press **Mode Setting** on the "setting item selection screen" (it returns to the preset tare input screen).

# Various settings (Setting methods (continued))

Eng

Before use  
(Various settings)


## Select the measurement flow (Continued from page 12).

**3** Press **8** and then press **Enter / Next**.

- The "ON or OFF of the one step mode" selection screen is displayed.

**Note**

- The one step mode is
  - a mode to measure body weight after inputting personal information.

**4** Set ON or OFF of the one step mode when using the body composition analyzer.


**Note**

- The default is "0. off." ("1. ON" for valid, "0. off" for invalid.)
- To correct the input value, or cancel the input,
  - press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.


- It returns to the "setting item selection" screen.

## Set ON or OFF of the athletic mode selection (Continued from page 12).

**3** Press **9** and then press **Enter / Next**.


- The "Set ON / OFF of the athletic mode selection" setting screen is displayed.

**4** Set ON or OFF of the athletic mode selection when using the body composition analyzer.


**Note**

- The default is "1. on." ("1. on" for valid, "0. off" for invalid.)
- To correct the input value, or cancel the input,
  - press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.

**Note**

- When various settings are continuously carried out,
  - press each number to set.
- The set contents are memorized until they are changed next time.


Eng

Before use  
(Various settings)

## Set the height increment (Continued from page 12).

**3** Press **1 0** and then press **Enter / Next**.


- The "height increment" selection screen is displayed.
- \* Invalid in "lb" mode.

**4** Selects the height increment.


**Note**

- The default is "1. on" (1 cm / 0.5 in).
  - 0. off: sets 0.1 cm increment
  - 1. on: sets 1 cm / 0.5 in increment
- To correct the input value, or cancel the input,
  - press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.


## Set the automatic determination time when inputting (Continued from page 12).

**3** Press **1 1** and then press **Enter / Next**.

- The "automatic determination time when inputting" setting screen is displayed.

**Note**


- Automatic determination when inputting is,
  - a function to determine the input value automatically even without pressing **Enter / Next**, after inputting the numeric value.

**4** Set the automatic determination time when inputting.


**Note**

- The default is 5 seconds ("5"). (input range: 0 - 9).
- \*If "0" is set, it is not automatically determined.
- To correct the input value, or cancel the input,
  - press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.


\*When various settings are all completed,
 

- press **Mode Setting** on the "setting item selection screen" (it returns to the preset tare input screen).

## Set ON / OFF of Wrestler mode (Continued from page 12).

**3** Press **1** **7** and then press **Enter / Next**.

**4** Set ON or OFF of Wrestler mode.

**Note**

- The default is 0.off
- 0: Set Wrestler mode off
- 5: Set Wrestler mode 5%
- 7: Set Wrestler mode 7%
- When the "Wrestler Mode" is activated, the SC-331S automatically calculates the Minimum Wrestling Weight at a predetermined minimum body fat %. The 1996 ACSM Position Stand "Weight Loss in Wrestlers" as adopted by the NCAA recommends the following MINIMUM body fat % as follows:


5% for Collegiate Athletes (male)

7% for High School Athletes (male)

**!** The target body fat % value is completely separate from the Minimum Wrestling Weight (Min Weight) calculations.


For example:

The target body fat % value may be entered as 15%, even though the Min Weight is calculated at a predetermined minimum body fat of 5 or 7%.

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.


**Note**

- \* When various settings are continuously carried out,  $\Rightarrow$ press each number to set.
- \* The set contents are memorized until they are changed next time.

## Set ON / OFF of the target body fat % (Continued from page 12).


**3** Press **1** **8** and then press **Enter / Next**.

**4** Set ON or OFF of the target body fat %.


**Note**

- The default is "0. off" ("1. on" for valid, "0. off" for invalid.)
- To correct the input value, or cancel the input,  $\Rightarrow$ press **CE** (the input is deleted).


**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.

## Set print language (Continued from page 12).


**3** Press **1** **9** and then press **Enter / Next**.

**4** Set the print language.


**Note**

- 1: English / 2: French / 5: Spanish / 7: Portuguese
- To correct the input value, or cancel the input,  $\Rightarrow$ press **CE** (the input is deleted).

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the "setting item selection" screen.

- \* When various settings are all completed,  $\Rightarrow$ press **Set Up** on the "setting item selection screen" (it returns to the preset tare input screen).

# Various settings Setting methods (continued)


Eng

Before use  
(Various settings)

## Set the print item preset (Continued from page 12).

**3** Press **2** **0** and then press **Enter / Next**.

- The “print out preset” setting screen is displayed.

**4** Set the print item preset.


- Set with **1** - **3**.

- Note**
- The print items preset are  $\Rightarrow$  a function to set the print items of the pattern 1 - 3 to ON automatically at once (page 23). The default is “1” (pattern 1).
  - When ON or OFF of print items are changed with the “setting of the print out item” after the print out preset (page 26), the state set for the last time is valid.

**5** After inputting the numeric value, press **Enter / Next**.

- It returns to the “setting item selection” screen.

**Note**

- \* When various settings are continuously carried out,  $\Rightarrow$  press each number to set.
- \* The set contents are memorized until they are changed next time.

## [Lists of contents of the print item preset]

Print item	Body composition analyzer									Scale
	1 Pattern 1			2 Pattern 2			3 Pattern 3			
	Standard	Athletic	Child	Standard	Athletic	Child	Standard	Athletic	Child	
Body type	✓	✓	✓	✓	✓	✓	✓	✓	✓	
TANITA Logo	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Category name	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Date and time	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Serial number	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Memo Space	✓	✓	✓	✓	✓	✓	✓	✓	✓	
ID No	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Fat mass	✓	✓	✓	✓	✓	✓				
Fat free mass	✓	✓	✓	✓	✓	✓				
Muscle mass	✓	✓		✓	✓					
Total body water	✓	✓		✓	✓					
Total body water %	✓	✓		✓	✓					
Bone mass	✓	✓		✓	✓					
BMR	✓	✓		✓	✓					
Metabolic age	✓	✓		✓	✓					
Visceral fat rating	✓	✓		✓	✓					
BMI	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Ideal body weight	✓			✓						
Degree of obesity	✓			✓						
Desirable range	✓	✓	✓							
Wrestler Mode										
Indicator Fat %	✓	✓	✓	✓	✓	✓				
Indicator BMI	✓	✓		✓	✓					
Indicator Visceral fat rating	✓	✓								
Indicator Muscle mass	✓	✓								
Indicator BMR	✓	✓								
Physique rating	✓	✓								

- The items marked with “✓” are printed.
- The items marked with “✓” can be selected to print (page 24)
- \* See page 24 for an example of preset print.

- \* When various settings are all completed,  $\Rightarrow$  press **Mode Setting** on the “setting item selection screen” (it returns to the preset tare input screen).

Eng

Before use  
(Various settings)

# Various settings Setting methods (continued)

Eng  
Before use  
(Various settings)

[In the case to select the print item preset "1"]

**TANITA**  
BODY COMPOSITION ANALYZER SC-331S

OCT/14/2009 18:21  
SERIAL No. 00000002

**LOGO**

**Serial No.**  
• The default is 00000001.  
Adds 1 each time it measures.

**ID**  
• When it is set with an ID, it is printed out. (The default is without an ID.)

**Fat %**  
• Fat % is amount of body fat as a proportion of body weight.

**FFM**  
• Fat Free Mass is comprised of muscle, bone, tissue, water, and all other fat free mass in the body.

**TBW\***  
• Total Body Water is the amount of water retained in the body. TBW is said to comprise between 50% - 70% of total body weight. Generally, men tend to have higher water weight than women due to a greater amount of muscle.

**Bone mass\***  
• Bone mineral amount included in the entire bone.

**Metabolic age\***  
• Metabolic age is evaluated young when a muscular amount is larger, and BMR is higher.

**BMI**  
• Calculated with "weight (kg) / height (m)<sup>2</sup>"  
• The standard value is for the Standard mode. In the case of the Athletic mode, the standard value is just a reference. And for those who are 17 years old or younger, only the body fat % is displayed as the standard value. The muscle mass, total body water and the estimated bone mass for those who are 17 years old or younger are for reference.

**Category name**

**Weight**  
• Measured weight.

**Fat mass**  
• Total weight of fat mass in the body.

**Muscle mass\***  
• Bone-free lean tissue mass (LTM)

**TBW %\***

**BMR\***  
• Basal Metabolic Rate represents the total energy expended by the body to maintain normal functions at rest such as respiration and circulation.

**Visceral fat rating\***  
• Visceral fat rating feature indicates the rating of visceral fat.

**Ideal body weight\***  
• Ideal body weight is a value for which the BMI is 22.

**Degree of obesity\***  
• Calculated as (weight - Ideal body weight) / Ideal body weight × 100.

**Impedance**  
• Impedance (This does not affect judgment of the measurement results.)

**INPUT**  
ID No. 000000123  
BODY TYPE STANDARD  
GENDER MALE  
AGE 36  
HEIGHT 170.4cm  
CLOTHES WEIGHT 1.0kg

**RESULT**  
WEIGHT 64.4kg  
FAT % 13.9 %  
FAT MASS 9.0kg  
FFM 55.4kg  
MUSCLE MASS 52.6kg  
TBW 39.2kg  
TBW % 60.9 %  
BONE MASS 2.8kg  
BMR 6678 kJ  
1596kcal  
METABOLIC AGE 21  
VISCERAL FAT RATING 4  
BMI 22.2  
IDEAL BODY WEIGHT 63.9kg  
DEGREE OF OBESITY 0.8 %

**DESIRED RANGE**  
FAT % 8.0-19.9 %  
FAT MASS 4.8-13.8kg

**Target**  
TARGET BF% is: 10 %  
Predicted weight: 61.6kg  
Predicted fat mass: 6.2kg  
FAT TO LOSE: 2.8kg

Consult your physician before beginning any weight management program. Tanita is not responsible for determining your targetBF%.

**INDICATOR**  
\*FAT %  
- | 0 | + | ++  
\*BMI  
- | 0 | + | ++  
\*VISCERAL FAT RATING  
| 13  
\*MUSCLE MASS  
- | 0 | +  
\*BMR  
- | 0 | +  
\*PHYSIQUE RATING  
STANDARD  
\*IMPEDANCE 481.2 Ω

\*18 - 99 years only

Please consult your doctor before you start a body weight management program. Tanita is not responsible for the target body fat ratio.

[In the case to select the print item preset "2"]

**TANITA**  
BODY COMPOSITION ANALYZER SC-331S

OCT/14/2009 18:22  
SERIAL No. 00000003

**INPUT**  
ID No. 000000123  
BODY TYPE STANDARD  
GENDER MALE  
AGE 36  
HEIGHT 170.4cm  
CLOTHES WEIGHT 1.0kg

**RESULT**  
WEIGHT 64.4kg  
FAT % 13.7 %  
FAT MASS 8.8kg  
FFM 55.6kg  
MUSCLE MASS 52.8kg  
TBW 39.5kg  
TBW % 61.3 %  
BONE MASS 2.8kg  
BMR 6699 kJ  
1601kcal  
METABOLIC AGE 21  
VISCERAL FAT RATING 4  
BMI 22.2  
IDEAL BODY WEIGHT 63.9kg  
DEGREE OF OBESITY 0.8 %

**Target**  
TARGET BF% is: 10 %  
Predicted weight: 61.8kg  
Predicted fat mass: 6.2kg  
FAT TO LOSE: 2.6kg

Consult your physician before beginning any weight management program. Tanita is not responsible for determining your targetBF%.

**INDICATOR**  
\*FAT %  
- | 0 | + | ++  
\*BMI  
- | 0 | + | ++  
\*IMPEDANCE 473.0 Ω

[In the case to select the print item preset "3"]

**TANITA**  
BODY COMPOSITION ANALYZER SC-331S

OCT/14/2009 18:23  
SERIAL No. 00000004

**INPUT**  
ID No. 000000123  
BODY TYPE STANDARD  
GENDER MALE  
AGE 36  
HEIGHT 170.4cm  
CLOTHES WEIGHT 1.0kg

**RESULT**  
WEIGHT 64.4kg  
FAT % 13.8 %  
BMI 22.2

**Target**  
TARGET BF% is: 10 %  
Predicted weight: 61.7kg  
Predicted fat mass: 6.2kg  
FAT TO LOSE: 2.7kg

Consult your physician before beginning any weight management program. Tanita is not responsible for determining your targetBF%.

\*IMPEDANCE 478.2 Ω

Please consult your doctor before you start a body weight management program. Tanita is not responsible for the target body fat ratio.

Eng  
Before use  
(Various settings)

# Various settings (Setting methods (continued))

Eng

Eng

Before use  
(Various settings)

Before use  
(Various settings)

## Set the items to print out (Continued from page 12).


**3** Select the number to set with the numeric keys ( **3** **9** - **6** **4** ) and press **Enter / Next** .

- The setting screen is displayed.

### Note

- See the next page for the setting number of each item.
- When the print item preset is set, after setting the items to print out (☞ page 24), the items to print may be changed automatically. Please confirm, "Print items preset list" (☞ page 23).

<In the case of selecting 39>

## Print item setting ON / OFF

39	TANITA Logo
40	Category name
41	Date and Time
42	Serial number
43	Memo space
44	ID Number
45	Fat mass
46	Fat free mass
47	Muscle mass
48	Total body water
49	Total body water %
50	Bone mass
51	Basal metabolic rate (BMR)


52	Metabolic age
53	Visceral fat rating
54	BMI
56	Ideal body weight
57	Degree of obesity
58	Desirable range Body fat %
59	Indicator - Fat %
60	Indicator - BMI
61	Indicator - Visceral fat rating
62	Indicator - Muscle mass
63	Indicator - BMR
64	Physique rating

**4** Selects ON or OFF of the item to print out.

### Note


- "0. off" is not to print out and "1. on" is to print out.
- To correct the input value, or cancel the input, ☞press **CE** (the input is deleted).

< In the case of selecting the logo setting on. >

**5** After inputting the numeric value, press **Enter / Next** .

- It returns to the "setting item selection" screen.

### Note

- \* When various settings are continuously carried out, ☞press each number to set.
- \* The set contents are memorized until they are changed next time.

- \* When various settings are all completed, ☞press **Mode Setting** on the "setting item selection screen" (it returns to the preset tare input screen).

**!** This explains the procedure when the printer is turned on. Please be aware that the display may be different if the number of sheets to print is set to **0**.

- !** Do not wipe the equipment with corrosive chemicals (gasoline, cleaner, etc.). Please use a neutral detergent to clean the equipment.
- When the equipment has been transferred to any location where there is a temperature difference of 36°F / 20°C or more, wait for at least two hours before using it.
- In taking measurements, please keep the person away from the unit, who uses transmitters such as a cellular phone avoid causing margin errors.
- The SC-331S estimates the various outputs specified in the Specification Section of this Manual. Keep in mind that there are many factors that can influence the accuracy of these measurements. Please refer to the General Instructions Section of this manual as this will help you minimize the variability (inaccuracy) of the measurements.

**About Athletic Mode**

- It is recommended for those who are 18 years old or older and meet the following conditions to select "Athletic Mode" and measure as reference values.
  - Those who exercise for 12 hours or more per a week.
  - Those who belong to a sport team or a sport organization with the aim of participation in competition, etc.
  - Those who exercise to build up like a bodybuilder.
  - Those who are professional athletes.

**Attention**

- The posture when measuring
  - Stand with both feet parallel on the electrodes.
  - Stand fullface without bending knees.
- The age input range is 5 – 99 years old.
  - Input age 99 for those who are 100 years old or older.

**Note**


- False results may be reported after excessive food/fluid intake, or after periods of intense exercise.
- When the user is 18 years or older, Athletic mode can be used.
- When Clothes weight is input, Clothes weight is subtracted from measurements and it displays it as Weight.
- "One step mode" should be the default.
- "Option mode" should move to back of "One step mode" page.

**In the case of the "one step mode"**

In the one step mode, after inputting the personal data, the weight and body composition are measured.

**1** Press **ON/OFF** to turn on the power.

(kg mode)

**2** Check that the body composition is selected, **input weight of clothing.** Input it by pressing **0** - **9** and **.**

**Note**



- The clothes weight (preset tare) can be input in the range of 0.0-20.0lb or 0.0-10.0kg.
  - \*0.2lb / 0.1kg increments
  - \*Unit is changed by **kg / lb** key.
- If it is mistakenly input,
  - >press **CE** (the input is deleted).

**3** Press **Enter / Next**.

**Note**


- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17). (The "body type selection" screen is displayed.)

**4** Input ID number. Input it by pressing **0** - **9**.

**Note**

- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17).
- The ID number can be input from 0 – 999999999. If **Enter / Next** is pressed, the non-inputted digits are filled with 0s.
- If it is mistakenly input,
  - >press **CE** (the input is deleted).
- If **CE** is pressed in the state that an ID number is not input, it returns to the "measurement start" screen.


## 5 Select body type.

Press the body type selection keys to input. When the body type is selected, the lamp flashes on the "Gender".

### Note

- This screen is not displayed if OFF is set in the "setting ON or OFF of the athletic mode selection" (page 18).
- The body type can also be selected with the numeric keys ( 1 2 ).
- If it is mistakenly input,
  - ⇨ press **CE** (the input is deleted, and it returns to the "body type selection" screen).
- If **CE** is pressed in the state that the body type is not input, it returns to the "ID number input" screen (or "measurement start" screen).


## 6 Select gender.

Press the male / female selection keys to input. When male or female is selected, the lamp flashes on the "Age".

### Note

- If it is mistakenly input,
  - ⇨ press **CE** (the input is deleted, and it returns to the "gender selection" screen).
- If **CE** is pressed in the state that the male or female is not selected, it returns to the "body type selection" screen.


## 7 Input age.

Input it by pressing 0 - 9.


### Note

- The age can be input from 5 - 99.
- If it is mistakenly input,
  - ⇨ press **CE** (the input is deleted).
- If **CE** is pressed in the state that the age is not input, it returns to the "gender selection" screen.

## 8 Press Enter / Next.

When the age is input, the lamp flashes on the "Height".

## 9 Input height.

Input it by pressing 0 - 9.

### Note

- The height can be input from 3ft -7ft11.5in / 90.0-249.9cm.
  - \*Unit is changed by **kg / lb** key.
- If it is mistakenly input,
  - ⇨ press **CE** (the input is deleted).
- If **CE** is pressed in the state that the height is not input, it returns to the "age input" screen.


## 10 Press Enter / Next.

The "Step on" lamp flashes.

### Note

- When **CE** is pressed, it returns to the previous screen.

## 11 Step on the electrodes with bare feet.

Take off your socks and stockings before stepping on.

## 12 Under measurement of the body composition.

88888 display goes off sequentially.

## 13 Measurement completion

The measurement result and the body fat percentage evaluation are displayed. It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (page 15).

### Note



- Fat % Healthy range is not displayed if OFF is set in the "setting ON or OFF of the Fat % Healthy range display" (page 17).

Step off the platform.

It returns to the "ID input" screen.

### Note

- The "body type selection" screen is displayed in the case that OFF is set in the "setting with or without an ID" (page 17).

**How to read the display <kg mode>**

Judgment result based on the body fat percentage. (page 42.)

**How to read the display <lb mode>**

Judgment result based on the body fat percentage. (page 42.)

## In the case of the "Option mode" (not the one step mode)


In the option mode, after measuring weight, personal data is input and then body composition is measured.

## 1 Press ON/OFF to turn on the power.

### Note

- Unit is changed by kg / lb key.


(kg mode)

## 2 Check that the body composition is selected and input weight of clothing.




Input it by pressing 0 - 9 and . .

### Note

- The clothes weight (preset tare) can be input in the range of 0.0-20.0lb or 0.0-10.0kg.  
\*0.2lb / 0.1kg increments
- To correct the input value,  
⇨press CE (the input is deleted).

## 3 Press Enter / Next .

The "Step on" lamp flashes.


### Note

- When CE is pressed, it returns to the previous screen.

## 4 Step on the electrodes with bare feet.


Take off your socks and stockings before stepping on.

When the weight becomes stable, it changes to the screen at the right.

### Note

- Do not step off the platform.
- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17). (The "body type selection" screen is displayed.)


## 5 Input ID number.

Input it by pressing 0 - 9.

### Note


- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17).
- The ID number can be input from 0 - 999999999. If **Enter / Next** is pressed, the non-inputted digits are filled with 0s.
- If it is mistakenly input, press **CE** (the input is deleted).
- If **CE** is pressed in the state that an ID number is not input, it returns to the "measurement start" screen.

## 6 Press **Enter / Next**.

### Note

- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17).


## 7 Select body type.

Press the body type selection keys to input. When the body type is selected, the lamp flashes on the "Gender".

### Note

- This screen is not displayed if OFF is set in the "setting ON or OFF of the athletic mode selection" (page 18).
- The body type can also be selected with the numeric keys (1 2).
- If it is mistakenly input, press **CE** (the input is deleted, and it returns to the "body type selection" screen).
- If **CE** is pressed in the state that the body type is not input, it returns to the "ID number input" screen (or "measurement start" screen).


## 8 Select gender.

Press the male / female selection keys to input. When male or female is selected, the lamp flashes on the "Age".

### Note

- If it is mistakenly input, press **CE** (the input is deleted, and it returns to the "gender selection" screen).
- If **CE** is pressed in the state that the male or female is not selected, it returns to the "body type selection" screen.


## 9 Input age.

Input it by pressing 0 - 9.

### Note


- The age can be input from 5 - 99.
- If it is mistakenly input, press **CE** (the input is deleted).
- If **CE** is pressed in the state that the age is not input, it returns to the "gender selection" screen.

## 10 Press **Enter / Next**.

When the age is input, the lamp flashes on the "Height".

## 11 Input height.

Input it by pressing **0** - **9**.

### Note

- The height can be input from 3ft -7ft11.5in / 90.0-249.9cm.
- \*Unit is changed by **kg / lb** key.
- If it is mistakenly input, press **CE** (the input is deleted).
- If **CE** is pressed in the state that the height is not input, it returns to the "age input" screen.

## 12 Press **Enter / Next**.

## 13 Under measurement of the body composition.

**88888** display goes off sequentially.

## 14 Measurement completion

The measurement result and the body fat percentage evaluation are displayed.



It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (page 15).

### Note

- Fat % Healthy range is not displayed if OFF is set in the "setting ON or OFF of the Fat % Healthy range display" (page 17).

Step off the platform.

It returns to the "measurement start" screen.

**How to read the display <kg mode>**

Judgment result based on the body fat percentage. (page 42.)

**How to read the display <lb mode>**

Judgment result based on the body fat percentage. (page 42.)

### 1 Input the target body fat %.

Input it by pressing **0** - **9**.

#### Note

- This screen is displayed if ON is set in the "setting ON or OFF of the target body fat %" (page 21).
- The target body fat % can be input from 4-55.
- If it is mistakenly input, press **CE** (the input is deleted).

If the number of sheet to print is set to "0", target body fat % setting function will be OFF, automatically.

If the target body fat is set to 0 or nothing, the target body fat % will not print.

**Caution** Before you start a body weight management program and set the appropriate personal body fat %, please consult your doctor. Tanita is not responsible for setting the appropriate target body fat % for specific individuals.

### 1 Press **ON/OFF** to turn on the power.

#### Note

- There is no automatic weight lock function.


### 2 Check that the body weight is selected and input the clothes weight.

Input it by pressing **0** - **9** and **.**

#### Note



- The clothes weight (preset tare) can be input in the range of 0.0-20.0lb or 0.0-10.0kg.  
\*0.2lb / 0.1kg increments
- \*Unit is changed by **kg / lb** key.
- If it is mistakenly input, press **CE** (the input is deleted).

### 3 Press **Enter / Next**.

#### Note

- This screen is not displayed if OFF is set in the "setting with or without an ID" (page 17). (The "measurement start" screen is displayed.)


## 4 Input an ID number.

Input it by pressing **0** - **9**.

### Note




- This screen is not displayed in the case that OFF is set in the "setting with or without an ID" (page 17).
- The ID number can be input from 0 - 9999999999.
- If it is mistakenly input, press **CE** (the input is deleted).
- If **CE** is pressed in the state that an ID number is not input, it returns to the "measurement start."

## 5 Press **Enter / Next**.

### Note

- This screen is not displayed in the case that OFF is set in the "setting with or without an ID" (page 17).

The lamp flashes for "Step on."


### Note

- If **CE** is pressed, it returns to the previous screen.

## 6 Measurement completion

It is automatically printed out. (In the case that other than 0 is set in the "Setting the number of sheets to print" (page 15).

Step off the platform.  
It returns to the "measurement start" screen.

## - What is body fat percentage? (Applicable age 5 - 99)

Body fat percentage is the amount of body fat as a proportion of your body weight. Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer. The chart below shows the healthy ranges for body fat.

Body Fat Ranges for Standard Children<sup>1</sup> <sup>1</sup>Susan Jebb et al. Obesity Research(02) 2004;12:A156-157  
 Body Fat Ranges for Standard Adults<sup>2</sup> <sup>2</sup>Gallagher D et al. Am J Clin Nutr 2000;72:694-701.  
<sup>1</sup>"New Body Fat Reference Curves for children"  
<sup>2</sup>"Healthy percentage body fat ranges: an approach for developing guidelines based on body mass index."

Age	Underfat				Healthy												Overfat				Obese																								
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
5	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
6	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
7	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
8	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
9	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
10	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
11	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
12	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
13	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
14	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
15	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
16	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
18	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
19	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
20-39	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
40-59	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
60-	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45

Age	0%				10%												20%				30%				40%																				
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19 <th>20</th> <th>21</th> <th>22</th> <th>23</th> <th>24</th> <th>25</th> <th>26</th> <th>27</th> <th>28</th> <th>29 <th>30</th><th>31</th><th>32</th><th>33</th><th>34</th><th>35</th><th>36</th><th>37</th><th>38</th><th>39 <th>40</th><th>41</th><th>42</th><th>43</th><th>44</th><th>45 </th></th></th>	20	21	22	23	24	25	26	27	28	29 <th>30</th> <th>31</th> <th>32</th> <th>33</th> <th>34</th> <th>35</th> <th>36</th> <th>37</th> <th>38</th> <th>39 <th>40</th><th>41</th><th>42</th><th>43</th><th>44</th><th>45 </th></th>	30	31	32	33	34	35	36	37	38	39 <th>40</th> <th>41</th> <th>42</th> <th>43</th> <th>44</th> <th>45 </th>	40	41	42	43	44	45
5	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
6	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
7	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
8	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
9	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
10	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
11	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
12	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
13	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
14	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
15	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
16	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
18	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
19	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45
20-39	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39</						

### - What is basal metabolic rate (BMR)? (Applicable age 18 - 99)

#### WHAT IS BMR?

Your Basal Metabolic Rate(BMR) is the minimum level of energy your body needs when at rest to function effectively including your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You burn calories when sleeping.

About 70% of calories consumed every day are used for your basal metabolism. In addition, energy is used when doing any kind of activity however; the more vigorous the activity is the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the quantity of muscles you have, therefor increasing your muscle mass will help increase your basal metabolism.

By studying healthy individuals, scientists have found that as people age, their metabolic rate changes. Basal metabolism rises as a child matures. After a peak at the age of 16 or 17, it typically starts to decrease gradually.

Having a higher basal metabolism will increase the number of calories used and help to decrease the amount of body fat. A low basal metabolic rate will make it harder to lose body fat and overall weight.

#### HOW DOES A TANITA BODY COMPOSITION ANALYZER CALCULATE BMR?

The basic way of calculating Basal Metabolic Rate BMR is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).\*

\* Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week : A Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition

### - What is metabolic age? (Applicable age 18 - 99)

This feature calculates your BMR and indicates the average age associated with that type of metabolism.

If your BMR Age is higher than your actual age, it is an indication that you need to improve your metabolic rate. Increased exercise will build healthy muscle tissue, which will improve your metabolic age.

You will obtain a reading between 12 and 90. Under 12 will be displayed as "12" and over 90 displayed as "90".

### - What is muscle mass? (Applicable age 18 - 99)

This feature indicates the weight of muscle in your body. the muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increase, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.


#### What is physique rating?

This feature assesses your physique according to the ratio of body fat ad muscle mass in your body. As you become more active and reduce the amount of body fat, your physique rating will also change accordingly. Even though your weight may not change, your muscle mass and body fat levels may be changing making you healthier and at lower risk of certain diseases.

Each person should set their own goal of which physique they would like and follow a diet and fitness program to meet that goal.

Result	Physique Rating	Explanation
1	Hidden obese	<b>Small Frame Obese</b>
		This person seems to have a healthy body type based on physical appearance; however, they have a high body fat % with low muscle mass level.
2	Obese	<b>Medium Frame Obese</b>
		This person has a high body fat percentage, with a moderate muscle mass level.
3	Solidly-built	<b>Large Frame Obese</b>
		This person has both a high body fat % and a high muscle mass.
4	Under exercised	<b>Low Muscle &amp; Average Body Fat%</b>
		This person has an average body fat % and a less than average muscle mass level
5	Standard	<b>Ave. Muscle &amp; Ave. Body Fat %</b>
		This person has average levels of both body fat and muscle mass.
6	Standard Muscular	<b>High Muscle &amp; Ave. Body Fat % (Athlete)</b>
		This person has an average body fat % and higher than normal muscle mass level.
7	Thin	<b>Low Muscle &amp; Low Fat</b>
		This person has both a lower than normal body fat % and muscle mass level.
8	Thin and muscular	<b>Thin and muscular (Athlete)</b>
		This person has lower than normal body fat % while having adequate muscle mass.
9	Very Muscular	<b>Very Muscular (Athlete)</b>
		This person has lower than normal body fat % while having above average muscle mass.

Source : Data from Columbia University (New York) & Tanita Institute (Tokyo)

## - What is bone mass? (Applicable age 18 - 99)

This feature indicates the amount of bone (bone mineral level, calcium or other minerals) in the body. Research has shown that exercise and the development of muscle tissue are related to stronger, healthier bones. While bone structure is unlikely to make noticeable changes in a short period, it is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. People worried about bone disease should consult their physician. People who suffer from osteoporosis or low bone densities due to advanced age, young age, pregnancy, hormonal treatment or other causes, may not get accurate estimations of their bone mass.

Below is the result of estimated bone masses of persons aged 20 to 40, who are said to have the largest amounts of bone masses, by weight. (Source : Tanita Body Weight Science Institute)

Please use the below charts as a guide to compare your bone mass reading.

### Women : Average of estimated bone mass

Weight (lb)			Weight (kg)		
Less than 110 lb	110lb - 165 lb	165 lb and up	Less than 50 kg	50 kg - 75 kg	75 kg and up
4.3 lb	5.3 lb	6.5 lb	1,95 kg	2,40 kg	2,95 kg

### Men : Average of estimated bone mass

Weight (lb)			Weight (kg)		
Less than 143 lb	143 lb - 209 lb	209 lb and up	Less than 65 kg	65 kg - 95 kg	95 kg and up
5.9 lb	7.3 lb	8.1 lb	2,66 kg	3,29 kg	3,69 kg

**Note:** Persons described below may obtain varying readings and should take the values given for reference purposes only.

- Elderly persons - Women during or after menopause
- People receiving hormone therapy




"Estimated bone mass" is a value estimated statistically based on its correlation with the fat-free amount (tissues other than the fat). "Estimated bone mass" does not give a direct judgment on the hardness or strength of the bones or the risks of bone fractures. If you have concerns over your bones, you are recommended to consult a specialist physician.

# Troubleshooting

Eng


- Please check the following before asking for repair.

When necessary  
(Troubleshooting)

	Symptom	Please check
How to measure	<b>Impedance measurement error</b> 	<ul style="list-style-type: none"> <li>• Measure with bare feet.</li> <li>• When the soles of your feet are dry, drop water for about 0.5 ml with the attached dropper on the electrodes before measurement.</li> <li>• Check the input information.</li> </ul>
	<b>Zero point error</b> 	<ul style="list-style-type: none"> <li>• Turn off the power, and remove the items on the platform and turn on the power again, and then redo the measurement.</li> </ul>
	<b>The measured weight is not stable.</b>	<ul style="list-style-type: none"> <li>• Is it installed at a place with vibrations?</li> <li>• Is the platform inclined? ⇒ Keep the platform horizontal. (☞ pages 8 and 9)</li> <li>• Is anything caught in the gaps of the platform? ⇒ Remove anything caught in the gaps.</li> </ul>
Display part	<b>Nothing is displayed even after turning on the power.</b>	<ul style="list-style-type: none"> <li>• Check that the power supply is connected correctly.</li> </ul>
	 is displayed.	<ul style="list-style-type: none"> <li>• The weight to measure exceeds the measurement range.</li> </ul>

Eng

When necessary  
(Troubleshooting)

	Symptom	Please check
Printer	<b>Run out of printer paper</b> 	<ul style="list-style-type: none"> <li>• Printer paper is not supplied. ⇒ Supply paper.</li> <li>⇒ In the case that the printer is not used, press <b>CE</b> and redo the initial setting.</li> </ul>
	<b>Printer cover open</b> 	<ul style="list-style-type: none"> <li>• The printer cover is open. ⇒ Properly close it.</li> <li>⇒ Check that the printer paper is not slanted.</li> </ul>
	<b>Paper does not come out.</b>	<p><b>Check the settings.</b></p> <ul style="list-style-type: none"> <li>• Is 0 set for the number of sheets to print in the "determination of the number of sheets to print?" ⇒ Press 1 - 3. (☞ page 15)</li> <li>• The printer may be broken. ⇒ Contact the agent from which you have purchased the product.</li> </ul>
	<b>The paper comes out, however, nothing has been printed.</b>	<ul style="list-style-type: none"> <li>• Is the reverse side of the printer paper set? ⇒ Set the paper properly. (☞ page 11)</li> <li>• The printer may be broken. ⇒ Contact the agent from which you have purchased the product.</li> </ul>

### Note

- In the case that an error is displayed other than the above, turn off the power once, and then measure again.  
If the same error is displayed repeatedly, contact our customer service center.

# Connection with a personal computer

 The USB interface enables input and output from this equipment.

## - Specifications

Communication standards	USB
Communication method	Asynchronous Transmission
Signal speed	9600 baud
Data bit length	8 bits
Parity	None
Stop bit	1 bit
Flow control	None
Terminator	CR+LF

### Note

- USB connector (B - Type 4 - pin female) are located on the back of the control box.
- Please provide your own cable as necessary as none are included.  
USB cable : A - Type 4 - pin (male) - B - Type 4 - pin (male)
- You must install the necessary driver onto your PC, available download from <http://www.tanita.com>.

## - Transmit data

 **Caution** Transmission Data is output immediately after measurement regardless of the status of the receiving equipment (personal computer, etc.). Therefore the receiving equipment needs to be ready to accept the data before measuring.

### (1) Output example (in the case to kg unit)

#### Note

- The data are divided with commas (,) for each data.
- The terminator (end of the data) is CR (ASCII code 0DH), LF (ASCII code 0AH).
- The control data for the items 1 – 4 are for expansion. They are not used presently so the receiving side can ignore them.
- Because of its specifications, the scale will output an irrelevant signal approximately 0.1 second after data output. Make sure you ignore this signal and do not import it.

#### Entire body data

Model number	Serial number	ID number	Date (mm/dd/yyyy)	Time					
MO	"XXXXXX"	SN	"XXXXXXXX"	ID	"XXXXXXXX"	Dt	"mm/dd/yyyy"	TI	"hh:mm"
Body type	Gender	Age	Height	Preset tare					
Bt	0or2	GE	1or2	AG	XX	Hm	XXX.X	Pt	XX.X
Body weight	Fat %	Fat mass	Fat free mass	Muscle mass					
Wk	XXX.X	FW	XX.X	rW	XXX.X	MW	XXX.X	mW	XXX.X
Muscle score	Bone mass	Body water	Body water %	BMI					
sW	XX	bW	XXX.X	wW	XXX.X	ww	XX.X	MI	XXX.X
Ideal body weight	Degree of obesity	Visceral fat rating	BMR (kJ)	BMR (kcal)					
Sw	XXX.X	OV	XXX.X	IF	XX	rb	XXXXX	rB	XXXXX
BMR score	Metabolic age	Target body fat %	Predicted weight	Predicted fat mass					
rJ	XX	rA	XX	gF	XX	gW	XXX.X	gf	XXX.X
Fat to gain / lose	Impedance (50kHz)	Checksum							
gt	XXXXX.X	ZF	XXXX.X	CS	XX				

## (2) Output data items

## &lt;kg / cm mode&gt;

Item	Header	format	kg / cm mode	Order of output			
				Body composition			Scale
				Adult	Athlete	Child	
Control data	{0		2 bytes fixed length	1	1	1	1
Control data	~0		1 byte fixed length	2	2	2	2
Control data	~1		1 byte fixed length	3	3	3	
Control data	~2		1 byte fixed length	4	4	4	
Model number	MO	"XXXXXX"	8 bytes fixed length ("SC-331")	5	5	5	3
Serial number	SN	"XXXXXXXX"	10 bytes fixed length	6	6	6	4
ID number	ID	"XXXXXXXXXX"	12 bytes fixed length	7	7	7	5
Date	Dt	"mm / dd / yyyy"	12 bytes fixed length ("mm / dd / yyyy")	8	8	8	6
Time	TI	"hh : mm"	7 bytes fixed length	9	9	9	7
Body type	Bt	0 or 2	1 byte fixed length (0:Standard 2:Athletic)	10	10	10	
Gender	GE	1 or 2	1 byte fixed length (1:Male 2:Female)	11	11	11	
Age	AG	XX	1-2 bytes variable length	12	12	12	
Height	Hm	XXX.X	4-5 bytes variable length (cm)	13	13	13	
Preset tare	Pt	XX.X	3-4 bytes variable length (increments 0.1kg)	14	14	14	8
Body weight	Wk	XXX.X	3-5 bytes variable length (increments 0.1kg)	15	15	15	9
Fat %	FW	XX.X	3-4 bytes variable length (increments 0.1%)	16	16	16	
Fat mass	fW	XXX.X	3-5 bytes variable length (increments 0.1kg)	17	17	17	
FFM	MW	XXX.X	3-5 bytes variable length (increments 0.1kg)	18	18	18	
Muscle mass	mW	XXX.X	3-5 bytes variable length (increments 0.1kg)	19	19		
Muscle score	sW	XX	1-2 bytes variable length (1-24)	20	20		
Bone mass	bW	XXX.X	3-5 bytes variable length (increments 0.1kg)	21	21		
Body water	wW	XXX.X	3-5 bytes variable length (increments 0.1kg)	22	22		
Body water %	ww	XX.X	4 bytes fixed length (increments 0.1%)	23	23		
BMI	MI	XXX.X	3-5 bytes variable length (increments 0.1)	24	24	19	
Ideal body weight	Sw	XXX.X	3-5 bytes variable length (increments 0.1kg)	25			
Degree of obesity	OV	XXX.X	3-5 bytes variable length (increments 0.1%)	26			
Visceral fat rating	IF	XX	1-2 bytes variable length	27	25		
BMR (kJ)	rb	XXXXX	1-5 bytes variable length (increments 1kJ)	28	26		
BMR (kcal)	rB	XXXXX	1-5 bytes variable length (increments 1kcal)	29	27		
BMR score	rJ	XX	1-2 bytes variable length (1-24)	30	28		
Metabolic age	rA	XX	2 bytes fixed length	31	29		
Target body fat %	gF	XX	1-2 bytes variable length (4-55)	32	30	20	
Predicted weight	gW	XXX.X	3-5 bytes variable length (increments 0.1kg)	33	31	21	
Predicted fat mass	gf	XXX.X	3-5 bytes variable length (increments 0.1kg)	34	32	22	
Fat gain / lose	gt	XXXXX.X	3-7 bytes variable length (increments 0.1kg)	35	33	23	
Impedance (50kHz)	ZF	XXXX.X	5-6 bytes variable length (increments 0.1)	36	34	24	
Checksum	CS	XX	2 bytes fixed length	37	35	25	10

## Output sample

- Oct-09-2009, 19:59, Standard, Male, 25 years old, Height 174.5cm  
 {0,16,~0,1,~1,1,~2,1,MO,"SC-331",SN,"0000000002",ID,"000000112",Dt,"10/09/2009",TI,"19:59",Bt,0,GE,1,AG,25,  
 Hm,174.5,Pt,1.5,Wk,65.6,FW,20.3,fW,13.3,MW,52.3,mW,49.6,sW,0,bW,2.7,wW,33.6,ww,51.2,MI,22.7,  
 Sw,63.6,OV,-5.8,IF,10,rb,7134,rB,1705,rJ,10,rA,30,gF,4,gW,123.4,gf,123.4,gt,-123.4,ZF,490.9,CS,C7

## &lt;lb / ft-in mode&gt;

Item	Header	format	lb / ft-in mode	Order of output			
				Body composition			Scale
				Adult	Athlete	Child	
Control data	{0		2 bytes fixed length	1	1	1	1
Control data	~0		1 byte fixed length	2	2	2	2
Control data	~1		1 byte fixed length	3	3	3	
Control data	~2		1 byte fixed length	4	4	4	
Model number	MO	"XXXXXX"	8 bytes fixed length ("SC-331")	5	5	5	3
Serial number	SN	"XXXXXXXX"	10 bytes fixed length	6	6	6	4
ID number	ID	"XXXXXXXXXX"	12 bytes fixed length	7	7	7	5
Date	Dt	"mm / dd / yyyy"	12 bytes fixed length ("mm / dd / yyyy")	8	8	8	6
Time	TI	"hh : mm"	7 bytes fixed length	9	9	9	7
Body type	Bt	0 or 2	1 byte fixed length (0:Standard 2:Athletic)	10	10	10	
Gender	GE	1 or 2	1 byte fixed length (1:Male 2:Female)	11	11	11	
Age	AG	XX	1-2 bytes variable length	12	12	12	
Height	Hi	XX.X	4 bytes variable length (in)	13	13	13	
Preset tare	Pt	XX.X	3-4 bytes variable length (increments 0.2lb)	14	14	14	8
Body weight	Wp	XXX.X	3-5 bytes variable length (increments 0.2lb)	15	15	15	9
Fat %	FW	XX.X	3-4 bytes variable length (increments 0.1%)	16	16	16	
Fat mass	fW	XXX.X	3-5 bytes variable length (increments 0.2lb)	17	17	17	
FFM	MW	XXX.X	3-5 bytes variable length (increments 0.2lb)	18	18	18	
Muscle mass	mW	XXX.X	3-5 bytes variable length (increments 0.2lb)	19	19		
Muscle score	sW	XX	1-2 bytes variable length (1-24)	20	20		
Bone mass	bW	XXX.X	3-5 bytes variable length (increments 0.2lb)	21	21		
Body water	wW	XXX.X	3-5 bytes variable length (increments 0.2lb)	22	22		
Body water %	ww	XX.X	4 bytes fixed length (increments 0.1%)	23	23		
BMI	MI	XXX.X	3-5 bytes variable length (increments 0.1)	24	24	19	
Ideal body weight	Sw	XXX.X	3-5 bytes variable length (increments 0.2lb)	25			
Degree of obesity	OV	XXX.X	3-5 bytes variable length (increments 0.1%)	26			
Visceral fat rating	IF	XX	1-2 bytes variable length	27	25		
BMR (kJ)	rb	XXXXX	1-5 bytes variable length (increments 1kJ)	28	26		
BMR (kcal)	rB	XXXXX	1-5 bytes variable length (increments 1kcal)	29	27		
BMR score	rJ	XX	1-2 bytes variable length (1-24)	30	28		
Metabolic age	rA	XX	2 bytes fixed length	31	29		
Target body fat %	gF	XX	1-2 bytes variable length (4-55)	32	30	20	
Predicted weight	gW	XXX.X	3-5 bytes variable length (increments 0.2lb)	33	31	21	
Predicted fat mass	gf	XXX.X	3-5 bytes variable length (increments 0.2lb)	34	32	22	
Fat gain / lose	gt	XXXXX.X	3-7 bytes variable length (increments 0.2lb)	35	33	23	
Impedance (50kHz)	ZF	XXXX.X	5-6 bytes variable length (increments 0.1)	36	34	24	
Checksum	CS	XX	2 bytes fixed length	37	35	25	10

## Body composition measurement by the BIA method.

### Introduction

This equipment provides estimated values for each measured value of body fat percentage, fat mass, fatfree mass, muscle mass and bone mass by the DXA method, estimated value for the total body water measured value by the dilution method and estimated value for the visceral fat rating by MRI method using the Bioelectrical Impedance Analysis (BIA method).

For measurement, a mode must be selected based on body type.

- 1) Standard (for 5-99 years of age)
- 2) Athletic (for Athletic persons who exercise considerably more than non-athlete)

Making a distinction by body type in the measurement mode produces more reliable body composition measurements for athletic persons, whose body compositions differ from those of average persons.


### - Principles of body composition measurement

BIA is a means of measuring body composition – fat mass, predicted muscle mass, etc. – by measuring bioelectrical impedance in the body. Fat within the body allows almost no electricity to pass through, while electricity passes rather easily through water, much of which is found in muscles. The degree of difficulty with which electricity passes through a substance is known as the electrical resistance, and the percentage of fat and other body constituents can be inferred from measurements of this resistance.

The Tanita Body Composition Analyzer measures body composition using a constant current source with a high frequency current (50kHz, 90µA). The 8 electrodes are positioned so that electric current is supplied from the electrodes on the tips of the toes of both feet, and voltage is measured on the heel of both feet. The current flows into the upper limbs or lower limbs, depending on the body part(s) to be measured.

### - What is the DXA method?

DXA was originally designed to measure bone mineral content, but in the full-body scan mode the body fat percentage, fat mass, and fat free mass of individual body parts (arms, legs, trunk) can also be measured. The image below shows one example of body composition measurement results obtained by DXA.

Body composition measurement results obtained by DXA (Lunar Co., Ltd; DPX-L)

### - What is dilution method?

In the dilution method, a labeled substance for a known amount is given and the concentration in equilibrium diffusing evenly is measured to obtain the total amount of the solvent that dilutes the labeled substance.

To measure the total body water (TBW), deuterium oxide (D<sub>2</sub>O) is generally used as the labeled substance. Deuterium oxide uses the overall total body water as dilution space so the total body water can be obtained. To obtain the extracellular fluid amount, sodium bromide (NaBr) is used as a labeled substance. Bromine (Br) is said to not enter the inside of cells, and uses extracellular fluid as the dilution space.


### - What is the visceral fat?

Visceral Adipose Tissue (VAT) has been associated with increased risk of developing lifestyle-related diseases. Accordingly, knowing and periodically checking the estimated VAT accumulation serves as one factor of a number of factors in assessing the prevention of lifestyle related diseases.


Tanita has developed the technology for measuring the VAT accumulation risk through bioelectrical impedance analysis (BIA) in comparison with image analysis applied to magnetic resonance imaging (MRI), in addition to the established technology for measuring the percent of body fat. The VAT accumulation risk is calculated by estimating the VAT area by the BIA method on the basis of MRI image processing. This method has a higher correlation than the estimation of the VAT accumulation risk based on BMI or abdominal circumference (waist circumference), allowing estimation that corresponds more precisely to individuals.

\*The VAT area by MRI is calculated by carrying out an image processing of the cross section of the lumbar vertebra L4-L5 regions.

(Fig. 1 - Fig. 3: Research results by N. Y. Columbia University and Jikei University Published by the North American Association for the Study of Obesity [NAASO] in 2004.)

<Fig. 3> Relationship between VAT Area by MRI and Estimated VAT Area by Tanita's BIA

## - Factors giving errors in measurement

In the BIA method, impedance is measured and the body composition is calculated based on the value. It is known that impedance changes by the amount of the total body water that occupies about 60% of weight and the change in its distribution and temperature change. Therefore, for the purpose of research or for daily repeating of measurements, the measurement conditions must be kept constant. Measurement under the changing conditions of temperature and total body water distribution or blood flow volume of extremities due to exercising, taking a bath, etc., affects the measurement result since the electric resistance in the body also changes.

Therefore, it is recommended to measure under the following conditions for stable measurement.

- 1) 3 hours have passed after getting up and normal lifestyle activities are carried out during this period. (The impedance transits staying at a high level if you remain sitting after getting up or drive a car, etc.)
- 2) 3 hours or more have passed after eating. (For 2 – 3 hours after eating, the impedance has a tendency to decrease.)
- 3) 12 hours or more have passed after vigorous exercise for measurement. (The tendency toward changes in impedance is not stable depending on the type and rigorousness of the exercise.)
- 4) If possible urinate before taking measurement.
- 5) For repeated measurements, measure at the same hour as much as possible. (At the same time of measurement of weight, the measurements can be made more stable by measuring at the same time of the day)

Very stable measured values can be obtained by measuring under the above conditions.

And in the development of this equipment, the following 6 items were set as conditions for the regression equation.

- 1) Prohibition of alcohol intake for 12 hours before measurement
- 2) Prohibition of excessive exercise for 12 hours before measurement.
- 3) Prohibition of excessive eating and drinking the day before measurement
- 4) Prohibition of eating and drinking for 3 hours before measurement
- 5) Urination just before measurement
- 6) Avoid measurements during menstruation (for female)

## 2) Inter-day changes

The diagrams below offer examples of actual measurements made of inter-day changes. A study was done to determine the degree of change in the impedance between the feet during dehydration; the first two days represent a normal daily routine, while in the latter two days a state of dehydration was induced using a sauna.


No significant inter-day change was measured in body weight, impedance between the feet, or body fat percentage during the normal daily routine. During the dehydrated state, however, a drop in body weight of 1kg was noted, with the impedance between the feet rising approximately  $15\Omega$  on the first day of dehydration and  $30\text{--}35\Omega$  on the second day. As a result, body fat percentage was up by around 1% on the first day of dehydration and by 1.5-2% on the second day.

As mentioned earlier, impedance increases when body weight is reduced (such as by dehydration), and decreases when body weight is increased through excess consumption of food and drink. The inter-day change in impedance is thus inversely proportional to the change in body weight.

These inter-day changes stem from such causes as:

- 1) Temporary increases in body weight (total body water) through overeating and overdrinking
- 2) Dehydration due to heavy sweating during vigorous exercise
- 3) Dehydration due to alcohol consumption or the use of diuretics
- 4) Dehydration due to heavy sweating during saunas, etc.

Accordingly, it is recommended that instructions be provided to the subject to help eliminate these causes when accurate measurements are needed.


# The Regression Formula for Basal Metabolic Rate (BMR)

It has long been said among medical and nutritional specialists that "The Basal Metabolic Rate (BMR) is more determined by the Fat Free Mass (FFM) than by the body weight" (Persons of a given body weight with a higher FFM will have a higher BMR), and that from the aspect of evaluating the body composition, should be estimated from the FFM. In addition, in cases of simple estimation formulae which can calculate from the height, weight and age, without evaluating the body composition, there was a problem with excessively high evaluations being given to obese persons with large body weight, and conversely excessively small BMR evaluations given to muscular athletes, though these are not as many in number. Currently, the BMR estimation recursion formula developed by Tanita, the manufacturer of body composition analyzers, based on their research, works by multiple regressive analysis using this FFM, and has a higher degree of accuracy in the individual differences in body composition. In order to derive the BMR, resting respiratory metabolism (Resting Energy Expenditure: REE) was measured using a breath gas analysis device, and this estimation recursion formula was created based on this data.

**<Figure 1> The Relationship Between Resting Energy Expenditure (REE) According to Breath Gas Analysis and Weight, FFM**

(Presented at Nutrition Week, Held in San Diego in 2002)

As shown in Figure 1: the REE (BMR) has a stronger relationship to the FFM than to body weight, and a difference is visible between males and females in the distribution trends. We see that in principle that we should calculate from the FFM rather than by the old formula centred on the relationship with weight.


**<Figure 2> Comparison of BMR Values from the TANITA Multiple regression model and Breath Analysis**

(Presented at Nutrition Week Held in San Diego in 2002)

The current BMR retrogression formula is a formula which acts on the principle of using the FFM value from the results of body composition measurement according to the BIA. A good relationship is shown in the BMR value based on actual breath analysis REE or  $R=0.9$  ( $p<0.0001$ ). These results were presented at the First Annual Nutrition Week (American College of Nutrition, American Society for Clinical Nutrition, American Society for Parenteral and Enteral Nutrition, North American Association for the Study of Obesity) held in 2002 in San Diego.

**NOTE:** This model has been calibrated for those between ages of 18-84. Those individuals outside of this age range may not be obtain accurate readings.

# Specifications

Eng

When necessary  
(Specifications)

<b>Model number</b>		SC-331S		
<b>Power source</b>		AC adapter (included) Center Minus		
		Input Voltage: 100-240 VAC 50/60 Hz 1.5A Output Voltage: 7 VDC		
		Rated Current: 4 A No Load Input Voltage: 7VDC		
<b>Power Consumption</b>		28 W		
<b>Impedance Measurement</b>	<b>Measurement System</b>	Tetra polar Bioelectrical Impedance Analysis		
	<b>Measurement Frequency</b>	50 kHz		
	<b>Measurement Current</b>	90 μA		
	<b>Electrode Material</b>	Pressure Contact Stainless Steel Foot Pads		
	<b>Measurement Style</b>	Between Both Feet		
	<b>Measurement Range</b>	150 - 1200Ω		
	<b>Accuracy at first calibration</b>	± 2%		
<b>Weight Measurement</b>	<b>Measurement System</b>	Strain Gauge Load Cell		
	<b>Maximum Capacity / Minimum Graduation</b>	270 kg / 0.1 kg 600 lb / 0.2 lb		
	<b>Accuracy at first calibration</b>	± 0.2 kg / ± 0.4 lb		
<b>Input Items</b>	<b>Clothes Weight</b>	0 - 10 kg / 0.1 kg increments	0 - 20 lb / 0.2 lb increments	
	<b>Gender</b>	Male / Female		
	<b>Body Type</b>	Standard (5 - 99 years) / Athletic (18 - 99 years)		
	<b>Age</b>	5 - 99 years old / 1 year increments		
	<b>Height</b>	90.0 - 249.9 cm / 0.1 cm increments	3 ft - 7 ft 11.5 in / 0.5 in increments	
	<b>Target Body Fat %</b>	4 - 55 %		

Eng

When necessary  
(Specifications)

<b>Output Items</b>	<b>Display</b>	<b>Target Body Fat %</b>	4 - 55 %		
		<b>Weight</b>	0 - 270 kg / 0.1 kg increments	600lb / 0.2 lb increments	
<b>Gender</b>		Male / Female			
<b>Body Type</b>		Standard / Athletic			
<b>Age</b>		5 - 99 years old / 1 year increments			
<b>Height</b>		90.0 - 249.9 cm / 0.1 cm increments	3 ft - 7 ft 11.5 in / 0.5 in increments		
	<b>FAT%</b>	3 - 75% / 0.1% increments			
	<b>Print-out</b>	<b>Logo</b>	TANITA LOGO (240x64 dot)		
		<b>Model number</b>	SC-331S		
		<b>Date and Time</b>	2009 / 1 / 1 - 2099 / 12 / 31		
		<b>Serial Number</b>	00000001 - 99999999		
		<b>ID</b>	000000000 - 999999999		
		<b>Body Type</b>	Standard (5 - 99 years) / Athletic (18 - 99 years)		
		<b>Gender</b>	Male / Female		
		<b>Age</b>	5 - 99 years old / 1 year increments		
		<b>Height</b>	90.0 - 249.9 cm / 0.1 cm increments	3 ft - 7 ft 11.5 in / 0.5 in increments	
		<b>Clothes Weight</b>	0 - 10 kg / 0.1 kg increments	0 - 20 lb / 0.2 lb increments	
		<b>Weight</b>	0 - 270 kg / 0.1 kg increments	600lb / 0.2 lb increments	
		<b>FAT%</b>	3 - 75% / 0.1% increments		
		<b>Fat Mass</b>	0.1 kg increments	0.2 lb increments	
		<b>FFM</b>	0.1 kg increments	0.2 lb increments	
		<b>Muscle Mass</b>	0.1 kg increments	0.2 lb increments	
		<b>TBW</b>	0.1 kg increments	0.2 lb increments	
		<b>TBW%</b>	15 - 85% / 0.1% increments		
		<b>Predicted Fat Mass</b>	0.1 kg increments	0.2 lb increments	
		<b>BMR</b>	1 kJ increments / 1 kcal increments		
		<b>Metabolic Age</b>	1 years increments ( 12 - 90 years)		
		<b>Visceral FAT Rating</b>	1 Level increments ( 1-59 Level)		
		<b>BMI</b>	0.1 increments		
		<b>Predicted Weight</b>	0.1 kg increments	0.2 lb increments	
		<b>Desirable range</b>	3 - 75% / 0.1% increments		
	<b>Indicator - FAT %</b>				
	<b>Indicator - BMI</b>				
	<b>Indicator - Visceral FAT Level</b>				
	<b>Indicator - Muscle Mass</b>				
	<b>Indicator - BMR</b>				
	<b>Physique Rating</b>				
	<b>Impedance</b>	150 - 1200Ω			
	<b>Display</b>	3 Rows, 5 Digits LCD			
	<b>Output Data Interface</b>	USB			
	<b>Temperature Range of Usage</b>	32°F to 95°F / 0°C to 35°C			
	<b>Relative Humidity</b>	30 - 80% (without condensation)			
	<b>Weight of Equipment (Remote Display Version)</b>	6.8 kg	15.0 lb		
	<b>Weight of Equipment (Column Mounted Version)</b>	12.1kg	26.7 lb		
<b>Size</b>	<b>Platform</b>	372 × 375 × 101 mm (14.6 × 14.8 × 4.0 in)			
	<b>Height (Column Mounted Version)</b>	1024 mm (40.3in)			

USA and Canada

#### Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and Canadian ICES-003. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

#### Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

### Disposal

This equipment is electronic device. Please dispose of this equipment appropriately as not the general household waste but electronic equipment. Please follow a regional regulation when you dispose of this.

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**Target**  
TARGET BF% is: 10 %  
Predicted weight: 61.8kg  
Predicted fat mass: 6.2kg  
FAT TO LOSE: 2.6kg  
Consult your physician before beginning any weight management program. Tanita is not responsible for determining your targetBF%.

**INDICATOR**  
\*FAT %  
- | 0 | + | ++  
\*BMI  
- | 0 | + | ++  
\*IMPEDANCE 473.0 Ω

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**RESULT**  
WEIGHT 64.4kg  
FAT % 13.8 %  
BMI 22.2

**Target**  
TARGET BF% is: 10 %  
Predicted weight: 61.7kg  
Predicted fat mass: 6.2kg  
FAT TO LOSE: 2.7kg  
Consult your physician before beginning any weight management program. Tanita is not responsible for determining your targetBF%.

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#### TANITA Europe B.V.
#### TANITA Health Equipment H.K.LTD.
#### TANITA India Private Limited
#### TANITA (Shanghai) Trading Co., Ltd.