



An in-house software generates segmentation images, which shows automated segmentation, manual segmentation and the difference map between the two segmentations to facilitate image analysts in visual qualitative assessment.


1012733, auto: LVEDV = 156mL, RVEDV = 147mL, LVM = 90g




1012733, man: LVEDV = 154mL, RVEDV = 151mL, LVM = 91g




1012733, diff: LVEDV_diff = 2mL, RVEDV_diff = -4mL, LVM_diff = -1g



1012733, auto: LVESV = 64mL, RVESV = 60mL



1012733, man: LVESV = 56mL, RVESV = 57mL



1012733, diff: LVESV_diff = 8mL, RVESV_diff = 3mL

