

Supplementary Table.1 Primer sequences used for real-time PCR assays

	Forward primer	Reverse primer
<i>Arnt</i>	5'-TCTCCCTCCCAGATGATGAC-3'	5'-CAATGTTGTGTCGGGAGATG-3'
<i>Ahr</i>	5'-CTGGTTGTCACAGCAGATGCCT-3'	5'-CGGTCTTCTGTATGGATGAGCTC-3'
<i>Ccl3</i>	5'-ACTGCCTGCTGCTTCTCCTACA-3'	5'-ATGACACCTGGCTGGGAGCAAA-3'
<i>Ccl5</i>	5'-CCTGCTGCTTTGCCTACCTCTC-3'	5'-ACACACTTGGCGGTTCCCTTCGA-3'
<i>Cxcl1</i>	5'-GCACCCAAACCGAAGTCATAG-3'	5'-AGAAGCCAGCGTTCACCAGA-3'
<i>Cxcl2</i>	5'-GCCCAGACAGAAGTCATAGCC-3'	5'-CTCCTCCTTTCCAGGTCAGTTA-3'
<i>Ahrr</i>	5'-ACATACGCCGGTAGGAAGAGA-3'	5'-GGTCCAGCTCTGTATTGAGGC-3'
<i>Cyp1a1</i>	5'-CATCACAGACAGCCTCATTGAGC-3'	5'-CTCCACGAGATAGCAGTTGTGAC-3'
<i>Cyp1a2</i>	5'-CATCCCCCACAGCACAACGA-3'	5'-GGTAAGAAACCGCTCTGGGC-3'
<i>Cyp1b1</i>	5'-GCCACTATTACGGACATCTTCGG-3'	5'-ACAACCTGGTCCAACCTCAGCCT-3'
<i>Cox2</i>	5'-GCGACATACTCAAGCAGGAGCA-3'	5'-AGTGGTAACCGCTCAGGTGTTG-3'
<i>Hprt</i>	5'-CTGGTGAAAAGGACCTCTCGAAG-3'	5'-CCAGTTTCACTAATGACACAAACG-3'