

Additional file 2. Detection of anti-pp65 reactivity from pp65<sub>386-403</sub>, pp65<sub>422-439</sub>, SA-C3d and PBS immunized sera. Immunoblot assay with sera from pp65<sub>386-403</sub> (*n*=9), pp65<sub>422-439</sub> (*n*=9), SA-C3d (*n*=5) and PBS (*n*=5) immunized mice against full-length pp65 at 8 weeks post immunization. Sera were used for tests at 1:250 dilution. P: positive control, 5000x diluted Rabbit anti-His tag antibody, N: negative control, 100x diluted normal serum. Secondary antibody: HRP conjugated anti-mouse/rabbit/human IgG antibody.