

Hu 306 with gene-centring							
IHC subtype	#number	Intrinsic subtype					
		Luminal A	Luminal B	Basal-like	Her2-enriched	Normal breast-like	Unclassified
Luminal A(ER+/HER2-)	73	55	16			2	
Luminal B(ER+/HER2+)	23	3	13		6		1
Basal-like(ER-/HER2-)	45		2	38	3	2	
HER2-enriched(ER-/HER2+)	28		1	3	20	4	

Hu 306 with DWD							
IHC subtype	#number	Intrinsic subtype					
		Luminal A	Luminal B	Basal-like	Her2-enriched	Normal breast-like	Unclassified
Luminal A(ER+/HER2-)	73	53	18			2	
Luminal B(ER+/HER2+)	23	4	14		4		1
Basal-like(ER-/HER2-)	45		2	38	3	2	
HER2-enriched(ER-/HER2+)	28		1	3	18	6	

PAM50 with gene-centring							
IHC subtype	#number	Intrinsic subtype					
		Luminal A	Luminal B	Basal-like	Her2-enriched	Normal breast-like	Unclassified
Luminal A(ER+/HER2-)	73	51	20			2	
Luminal B(ER+/HER2+)	23	3	14		6		
Basal-like(ER-/HER2-)	45		2	37	4		2
HER2-enriched(ER-/HER2+)	28	2		4	20		2

PAM50 with DWD							
IHC subtype	#number	Intrinsic subtype					
		Luminal A	Luminal B	Basal-like	Her2-enriched	Normal breast-like	Unclassified
Luminal A(ER+/HER2-)	73	43	11			19	
Luminal B(ER+/HER2+)	23	6	6		11		
Basal-like(ER-/HER2-)	45	1	1	37	4		2
HER2-enriched(ER-/HER2+)	28	1	1	4	16		6