

Additional file 1: Neutralization of CCL3 by a specific antibody. Astrocytes primed with $1\mu g/ml$ LPS for 3h were treated with $10\mu M$ A β 42 for 3h. Specific antibody against CCL3 was added 45min after A β 42. The concentration of CCL3 in culture supernatants were then measured by ECLIA. Results are the mean \pm SEM derived from 3 experiments in duplicate and statistical differences were determined using a Mann Whitney test. ***P<0.001, **P<0.01 and *P<0.05 vs respective controls.