

Active sensitization

Cyclophosphamide
30mg/kg ip



OVA + Al(OH)₃ ip.
2mg 100mg



OVA + Al(OH)₃ ip.
0.01mg 100mg



-2day

0

3W

Measurement of
cough reflex sensitivity

Procaterol ip.
0.1mg/kg



0.5h

OVA inhalation
challenge




0.5h



6W+1day

First vehicle
or ONO 1714

ENO measurement
and BAL



0.5h

6W+2days

Second vehicle
or ONO 1714

